

Perforin-1

Catalog No: #AP79809



Package Size: #AP79809-1 50ug #AP79809-2 100ug #AP79809-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Perforin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	P1,Cytolysin,Lymphocyte pore-forming protein,PFP,PFP
Accession No.	Uniprot ID:P14222
Uniprot	P14222
GeneID	5551;
SDS-PAGE MW	33.99kDa
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Expression Region:247-555aaStorage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PRF1

Note: This product is for in vitro research use only