

Adenosine deaminase domain-containing protein 1

Catalog No: #AP79812



Package Size: #AP79812-1 50ug #AP79812-2 100ug #AP79812-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Adenosine deaminase domain-containing protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Testis nuclear RNA-binding protein, Tenr
Accession No.	Uniprot ID:Q3KR54
Uniprot	Q3KR54
GeneID	294979;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Adad1

Note: This product is for in vitro research use only