

## Norrin

Catalog No: #AP79813



Package Size: #AP79813-1 50ug #AP79813-2 100ug #AP79813-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Norrin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Norrie disease protein,X-linked exudative vitreoretinopathy 2 protein,EVR2
Accession No.	Uniprot ID:Q00604
Uniprot	Q00604
GeneID	4693;
SDS-PAGE MW	11.99kDa
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Expression Region:25-133aaStorage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:NDP

Note: This product is for in vitro research use only