

Nuclear receptor coactivator 3

Catalog No: #AP79814



Package Size: #AP79814-1 50ug #AP79814-2 100ug #AP79814-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Nuclear receptor coactivator 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	NCoA-3,ACTR,Amplified in breast cancer 1 protein,AIB-1,CBP-interacting protein,pCIP,Class E basic helix-loop-helix protein 42,bHLHe42,Receptor-associated coactivator 3,RAC-3,Steroid receptor coactivator protein 3,SRC-3,Thyroid hormone receptor activator molecule 1,TRAM-1,AIB1,BHLHE42,RAC3,TRAM1
Accession No.	Uniprot ID:Q9Y6Q9
Uniprot	Q9Y6Q9
GeneID	8202;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:NCOA3

Note: This product is for in vitro research use only