

P2Y purinoceptor 12

Catalog No: #AP79815



Package Size: #AP79815-1 50ug #AP79815-2 100ug #AP79815-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	P2Y purinoceptor 12
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	P2Y12,P2Y12 platelet ADP receptor,P2y12
Accession No.	Uniprot ID:Q9EPX4
Uniprot	Q9EPX4
GeneID	64803;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:P2ry12

Note: This product is for in vitro research use only