

Proteasome subunit beta type-9

Catalog No: #AP79817



Package Size: #AP79817-1 50ug #AP79817-2 100ug #AP79817-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Proteasome subunit beta type-9
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Low molecular mass protein 2,Macropain chain 7,Multicatalytic endopeptidase complex chain 7,Proteasome chain 7,Proteasome subunit beta-1i,Really interesting new gene 12 protein,LMP2,PSMB6i,RING12
Accession No.	Uniprot ID:P28065
Uniprot	P28065
GeneID	5698;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PSMB9

Note: This product is for in vitro research use only