

6-phosphofructo-2-kinase/fructose-2,6-bisphosphatase 4

Catalog No: #AP79818

Package Size: #AP79818-1 50ug #AP79818-2 100ug #AP79818-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	6-phosphofructo-2-kinase/fructose-2,6-bisphosphatase 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	6PF-2-K/Fru-2,6-P2ase 4,6PF-2-K/Fru-2,6-P2ase testis-type isozyme
Accession No.	Uniprot ID:Q16877
Uniprot	Q16877
GeneID	5210;
SDS-PAGE MW	12.1kDa
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Expression Region:258-367aaStorage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PFKFB4

Note: This product is for in vitro research use only