

Serine/threonine-protein kinase PINK1, mitochondrial

Catalog No: #AP79819



Package Size: #AP79819-1 50ug #AP79819-2 100ug #AP79819-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Serine/threonine-protein kinase PINK1, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	BRPK,PTEN-induced putative kinase protein 1
Accession No.	Uniprot ID:Q9BXM7
Uniprot	Q9BXM7
GeneID	65018;
SDS-PAGE MW	22.66kDa
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Expression Region:260-465aaStorage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PINK1

Note: This product is for in vitro research use only