

Promotilin

Catalog No: #AP79820



Package Size: #AP79820-1 50ug #AP79820-2 100ug #AP79820-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Promotilin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Accession No.	Uniprot ID:P12872
Uniprot	P12872
GeneID	4295;
SDS-PAGE MW	9.9kDa
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Expression Region:26-115aaStorage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MLN

Note: This product is for in vitro research use only