

## N-fatty-acyl-amino acid synthase/hydrolase PM20D1

Catalog No: #AP79821



Package Size: #AP79821-1 50ug #AP79821-2 100ug #AP79821-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	N-fatty-acyl-amino acid synthase/hydrolase PM20D1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Peptidase M20 domain-containing protein 1
Accession No.	Uniprot ID:Q6GTS8
Uniprot	Q6GTS8
GeneID	148811;
SDS-PAGE MW	25.85kDa
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Expression Region:26-260aaStorage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:PM20D1

Note: This product is for in vitro research use only