

Angiotensin-converting enzyme 2

Catalog No: #AP79823

Package Size: #AP79823-1 50ug #AP79823-2 100ug #AP79823-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Angiotensin-converting enzyme 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Angiotensin-converting enzyme 2,Angiotensin-converting enzyme 2,Agf
Accession No.	Uniprot ID:Q8R0Z6
Uniprot	Q8R0Z6
GeneID	70726;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Angptl6

Note: This product is for in vitro research use only