

# Squalene monooxygenase

Catalog No: #AP79824



Package Size: #AP79824-1 50ug #AP79824-2 100ug #AP79824-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Squalene monooxygenase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Squalene epoxidase,SE,ERG1
Accession No.	Uniprot ID:Q14534
Uniprot	Q14534
GeneID	6713;
SDS-PAGE MW	21.23kDa
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Expression Region:268-460aaStorage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:SQLE

Note: This product is for in vitro research use only