

# V-set and transmembrane domain-containing protein 2-like protein

Catalog No: #AP79826

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP79826-1 50ug #AP79826-2 100ug #AP79826-3 1mg

## Description

Product Name	V-set and transmembrane domain-containing protein 2-like protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Gm691
Accession No.	Uniprot ID:Q6PDS0
Uniprot	Q6PDS0
GeneID	277432;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Vstm2l

Note: This product is for in vitro research use only