

Insulin-like growth factor-binding protein 7

Catalog No: #AP79828



Package Size: #AP79828-1 50ug #AP79828-2 100ug #AP79828-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Insulin-like growth factor-binding protein 7
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	IBP-7,IGFBP-rP1,MAC25 protein,PGI2-stimulating factor,Prostacyclin-stimulating factor,Tumor-derived adhesion factor,TAF,MAC25,PSF
Accession No.	Uniprot ID:Q16270
Uniprot	Q16270
GeneID	3490;
SDS-PAGE MW	28.16kDa
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Expression Region:27-282aaStorage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:IGFBP7

Note: This product is for in vitro research use only