

Beta-1,3-galactosyltransferase 2

Catalog No: #AP79830



Package Size: #AP79830-1 50ug #AP79830-2 100ug #AP79830-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Beta-1,3-galactosyltransferase 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Beta-1,3-GalTase 2,UDP-galactose:2-acetamido-2-deoxy-D-glucose 3beta-galactosyltransferase 2
Accession No.	Uniprot ID:O43825
Uniprot	O43825
GeneID	8707;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:B3GALT2

Note: This product is for in vitro research use only