

Bile acyl-CoA synthetase

Catalog No: #AP79833



Package Size: #AP79833-1 50ug #AP79833-2 100ug #AP79833-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bile acyl-CoA synthetase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	BACS,Bile acid-CoA ligase,BA-CoA ligase,Cholate--CoA ligase,Fatty acid transport protein 5,FATP-5,Solute carrier family 27 member 5,Very long-chain acyl-CoA synthetase-related protein,VLACS-related,Acsb,Acsvl6,Fatp5,Vlacsr
Accession No.	Uniprot ID:Q4LDG0
Uniprot	Q4LDG0
GeneID	26459;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Slc27a5

Note: This product is for in vitro research use only