

Interferon-induced transmembrane protein 10

Catalog No: #AP79834



Package Size: #AP79834-1 50ug #AP79834-2 100ug #AP79834-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Interferon-induced transmembrane protein 10
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Dispanin subfamily A member 3,DSPA3
Accession No.	Uniprot ID:Q8BR26
Uniprot	Q8BR26
GeneID	320802;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Ifitm10

Note: This product is for in vitro research use only