

Immunoglobulin J chain

Catalog No: #AP79835



Package Size: #AP79835-1 50ug #AP79835-2 100ug #AP79835-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Immunoglobulin J chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Joining chain of multimeric IgA and IgM,IGCJ,IGJ
Accession No.	Uniprot ID:P01591
Uniprot	P01591
GeneID	3512;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:JCHAIN

Note: This product is for in vitro research use only