

Osteopontin

Catalog No: #AP79836



Package Size: #AP79836-1 50ug #AP79836-2 100ug #AP79836-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Osteopontin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	2AR,Bone sialoprotein 1,Calcium oxalate crystal growth inhibitor protein,Early T-lymphocyte activation 1 protein,Minopontin,Secreted phosphoprotein 1,SPP-1,Eta-1,Op,Spp-1
Accession No.	Uniprot ID:P10923
Uniprot	P10923
GeneID	20750;
SDS-PAGE MW	10.45kDa
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Expression Region:29-123aaStorage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Spp1

Note: This product is for in vitro research use only