

## CD1 Monoclonal Antibody(9H6)

Catalog No: #44036

Package Size: #44036-1 50ul #44036-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	CD1 Monoclonal Antibody(9H6)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	9H6
Isotype	IgG1
Applications	IHC
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous CD1 proteins.
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide
Target Name	CD1
Accession No.	Swiss-Prot:P06126Gene ID:909
Uniprot	P06126
GeneID	909;
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

Immunohistochemistry: 1:50-1:200

## Images



Immunohistochemical analysis of paraffin-embedded human-tonsils using antibody diluted at 1:50.

Note: This product is for in vitro research use only