

## Protein FAM132B

Catalog No: #AP79840



Package Size: #AP79840-1 50ug #AP79840-2 100ug #AP79840-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Protein FAM132B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Complement C1q tumor necrosis factor-related protein 15, Myonectin, C1QTNF15, CTRP15
Accession No.	Uniprot ID:Q4G0M1
Uniprot	Q4G0M1
GeneID	151176;
SDS-PAGE MW	13.42kDa
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Expression Region:29-150aa Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name:FAM132B

Note: This product is for in vitro research use only