

## Growth/differentiation factor 11

Catalog No: #AP79841



Package Size: #AP79841-1 50ug #AP79841-2 100ug #AP79841-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Growth/differentiation factor 11
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	GDF-11,Bone morphogenetic protein 11,BMP-11,Bmp11
Accession No.	Uniprot ID:Q9Z217
Uniprot	Q9Z217
GeneID	29454;
SDS-PAGE MW	11.99kDa
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Expression Region:297-405aaStorage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Gdf11

Note: This product is for in vitro research use only