

Probable phosphoglycerate mutase 4

Catalog No: #AP79844

Package Size: #AP79844-1 50ug #AP79844-2 100ug #AP79844-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Probable phosphoglycerate mutase 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	PGAM3
Accession No.	Uniprot ID:Q8N0Y7
Uniprot	Q8N0Y7
GeneID	441531;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PGAM4

Note: This product is for in vitro research use only