

## Dysbindin

Catalog No: #AP79847



Package Size: #AP79847-1 50ug #AP79847-2 100ug #AP79847-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Dysbindin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Biogenesis of lysosome-related organelles complex 1 subunit 8,BLOC-1 subunit 8,Dysbindin-1,Dystrobrevin-binding protein 1,Hermansky-Pudlak syndrome 7 protein homolog,HPS7 protein homolog,Bloc1s8,Sdy
Accession No.	Uniprot ID:Q91WZ8
Uniprot	Q91WZ8
GeneID	94245;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Dtnbp1

Note: This product is for in vitro research use only