

Elongation of very long chain fatty acids protein 1

Catalog No: #AP79848

Package Size: #AP79848-1 50ug #AP79848-2 100ug #AP79848-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Elongation of very long chain fatty acids protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	3-keto acyl-CoA synthase Elov11,ELOVL fatty acid elongase 1,ELOVL FA elongase 1,Very long chain 3-ketoacyl-CoA synthase 1,Very long chain 3-oxoacyl-CoA synthase 1,Ssc1
Accession No.	Uniprot ID:Q9JLJ5
Uniprot	Q9JLJ5
GeneID	54325;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Elov11

Note: This product is for in vitro research use only