

# Carboxyl-terminal PDZ ligand of neuronal nitric oxide synthase protein



Catalog No: #AP79849

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP79849-1 50ug #AP79849-2 100ug #AP79849-3 1mg

## Description

Product Name	Carboxyl-terminal PDZ ligand of neuronal nitric oxide synthase protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	C-terminal PDZ ligand of neuronal nitric oxide synthase protein, Nitric oxide synthase 1 adaptor protein, Capon
Accession No.	Uniprot ID:O54960
Uniprot	O54960
GeneID	192363;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Nos1ap

Note: This product is for in vitro research use only