

Vitelline membrane outer layer protein 1 homolog

Catalog No: #AP79850



Package Size: #AP79850-1 50ug #AP79850-2 100ug #AP79850-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Vitelline membrane outer layer protein 1 homolog
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Gm741
Accession No.	Uniprot ID:Q5SXG7
Uniprot	Q5SXG7
GeneID	327956;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Vmo1

Note: This product is for in vitro research use only