

Agrin

Catalog No: #AP79851



Package Size: #AP79851-1 50ug #AP79851-2 100ug #AP79851-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Agrin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	AGRIN
Accession No.	Uniprot ID:O00468
Uniprot	O00468
GeneID	375790;
SDS-PAGE MW	16.39kDa
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Expression Region:311-459aaStorage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:AGRN

Note: This product is for in vitro research use only