

Serine protease 29

Catalog No: #AP79854



Package Size: #AP79854-1 50ug #AP79854-2 100ug #AP79854-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Serine protease 29
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Implantation serine proteinase 2,ISP-2,Strypsin-2,Strypsin-related protein,Tryptase-like proteinase,Isp2
Accession No.	Uniprot ID:Q99MS4
Uniprot	Q99MS4
GeneID	114662;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Prss29

Note: This product is for in vitro research use only