

Protein-L-isoaspartate(D-aspartate) O-methyltransferase

Catalog No: #AP79857

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP79857-1 50ug #AP79857-2 100ug #AP79857-3 1mg

Description

Product Name	Protein-L-isoaspartate(D-aspartate) O-methyltransferase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	PIMT,L-isoaspartyl protein carboxyl methyltransferase,Protein L-isoaspartyl/D-aspartyl methyltransferase,Protein-beta-aspartate methyltransferase
Accession No.	Uniprot ID:P15246
Uniprot	P15246
GeneID	613854;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PCMT1

Note: This product is for in vitro research use only