

# Pleckstrin homology domain-containing family O member 1

Catalog No: #AP79858

Package Size: #AP79858-1 50ug #AP79858-2 100ug #AP79858-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Pleckstrin homology domain-containing family O member 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	PH domain-containing family O member 1,RCJMB04_21c20
Accession No.	Uniprot ID:Q5F3C8
Uniprot	Q5F3C8
GeneID	426108;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:PLEKHO1

Note: This product is for in vitro research use only