

Pleckstrin homology domain-containing family O member 1

Catalog No: #AP79859

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP79859-1 50ug #AP79859-2 100ug #AP79859-3 1mg

Description

Product Name	Pleckstrin homology domain-containing family O member 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	PH domain-containing family O member 1,C-Jun-binding protein,JBP,Casein kinase 2-interacting protein 1,CK2-interacting protein 1,Osteoclast maturation-associated gene 120 protein,CKIP1,OC120,HQ0024c
Accession No.	Uniprot ID:Q53GL0
Uniprot	Q53GL0
GeneID	51177;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PLEKHO1

Note: This product is for in vitro research use only