

Tafazzin

Catalog No: #AP79860

Package Size: #AP79860-1 50ug #AP79860-2 100ug #AP79860-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tafazzin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Protein G4.5,EFE2,G4.5
Accession No.	Uniprot ID:Q16635
Uniprot	Q16635
GeneID	6901;
SDS-PAGE MW	28.82kDa
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Expression Region:31-292aaStorage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TAZ

Note: This product is for in vitro research use only