

Lysine-specific histone demethylase 1A

Catalog No: #AP79861



Package Size: #AP79861-1 50ug #AP79861-2 100ug #AP79861-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Lysine-specific histone demethylase 1A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	BRAF35-HDAC complex protein BHC110,Flavin-containing amine oxidase domain-containing protein 2,Aof2,Kiaa0601,Lsd1
Accession No.	Uniprot ID:Q6ZQ88
Uniprot	Q6ZQ88
GeneID	99982;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Kdm1a

Note: This product is for in vitro research use only