

Kelch-like ECH-associated protein 1

Catalog No: #AP79862



Package Size: #AP79862-1 50ug #AP79862-2 100ug #AP79862-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Kelch-like ECH-associated protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	Cytosolic inhibitor of Nrf2, INrf2, INRF2
Accession No.	Uniprot ID:Q684M4
Uniprot	Q684M4
GeneID	100136900;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: KEAP1

Note: This product is for in vitro research use only