

Solute carrier family 12 member 1

Catalog No: #AP79864



Package Size: #AP79864-1 50ug #AP79864-2 100ug #AP79864-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Solute carrier family 12 member 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Bumetanide-sensitive sodium-(potassium)-chloride cotransporter 2, Kidney-specific Na-K-Cl symporter, NKCC2
Accession No.	Uniprot ID:Q13621
Uniprot	Q13621
GeneID	6557;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: SLC12A1

Note: This product is for in vitro research use only