

## DNA mismatch repair protein Msh2

Catalog No: #AP79866



Package Size: #AP79866-1 50ug #AP79866-2 100ug #AP79866-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	DNA mismatch repair protein Msh2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	hMSH2,MutS protein homolog 2
Accession No.	Uniprot ID:P43246
Uniprot	P43246
GeneID	4436;
SDS-PAGE MW	27.94kDa
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Expression Region:321-574aaStorage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:MSH2

Note: This product is for in vitro research use only