

Histone deacetylase 2

Catalog No: #AP79868



Package Size: #AP79868-1 50ug #AP79868-2 100ug #AP79868-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Histone deacetylase 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	HD2
Accession No.	Uniprot ID:Q92769
Uniprot	Q92769
GeneID	3066;
SDS-PAGE MW	17.38kDa
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Expression Region:331-488aaStorage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:HDAC2

Note: This product is for in vitro research use only