

## CD5 Monoclonal Antibody(10G8)

Catalog No: #44039

Package Size: #44039-1 50ul #44039-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	CD5 Monoclonal Antibody(10G8)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	10G8
Isotype	IgG1
Applications	IHC
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous CD5 proteins.
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide
Target Name	CD5
Other Names	CD5; LEU1; T-cell surface glycoprotein CD5; Lymphocyte antigen T1/Leu-1; CD5
Accession No.	Swiss-Prot:P06127Gene ID:921
Uniprot	P06127
GeneID	921;
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

Immunohistochemistry: 1:50-1:200

## Images



Immunohistochemical analysis of paraffin-embedded Human Tonsil Carcinoma using CD5 Mouse mAb diluted at 1:200.

---

Note: This product is for in vitro research use only