

## Orexin receptor type 2

Catalog No: #AP79871



Package Size: #AP79871-1 50ug #AP79871-2 100ug #AP79871-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Orexin receptor type 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Ox-2-R,Ox2-R,Ox2R,Hypocretin receptor type 2
Accession No.	Uniprot ID:Q9TUP7
Uniprot	Q9TUP7
GeneID	399545;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:HCRTR2

Note: This product is for in vitro research use only