

## Leukemia inhibitory factor receptor

Catalog No: #AP79872



Package Size: #AP79872-1 50ug #AP79872-2 100ug #AP79872-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Leukemia inhibitory factor receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	LIF receptor,CD118
Accession No.	Uniprot ID:P42702
Uniprot	P42702
GeneID	3977;
SDS-PAGE MW	29.7kDa
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Expression Region:335-604aaStorage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:LIFR

Note: This product is for in vitro research use only