

Interleukin-27 subunit beta

Catalog No: #AP79874



Package Size: #AP79874-1 50ug #AP79874-2 100ug #AP79874-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Interleukin-27 subunit beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	IL-27 subunit beta, Epstein-Barr virus-induced gene 3 protein, EBV-induced gene 3 protein, IL27B
Accession No.	Uniprot ID:Q14213
Uniprot	Q14213
GeneID	10148;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: EB13

Note: This product is for in vitro research use only