

Kelch-like protein 22

Catalog No: #AP79875



Package Size: #AP79875-1 50ug #AP79875-2 100ug #AP79875-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Kelch-like protein 22 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Hu |
| Accession No. | Uniprot ID:Q53GT1 |
| Uniprot | Q53GT1 |
| GeneID | 84861; |
| SDS-PAGE MW | 23.1kDa |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Expression Region:371-580aaStorage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:KLHL22

Note: This product is for in vitro research use only