

Kallikrein-1

Catalog No: #AP79876



Package Size: #AP79876-1 50ug #AP79876-2 100ug #AP79876-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Kallikrein-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Glandular kallikrein K1, KAL-B, Renal kallikrein, Tissue kallikrein-6, mGK-6, Kik-6, Kik6
Accession No.	Uniprot ID:P15947
Uniprot	P15947
GeneID	16612;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Kik1

Note: This product is for in vitro research use only