

Amyloid beta A4 protein

Catalog No: #AP79877



Package Size: #AP79877-1 50ug #AP79877-2 100ug #AP79877-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Amyloid beta A4 protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	ABPP,APP,Alzheimer disease amyloid A4 protein homolog,Amyloid precursor protein,Amyloidogenic glycoprotein,AG,Beta-amyloid precursor protein
Accession No.	Uniprot ID:P12023
Uniprot	P12023
GeneID	11820;
SDS-PAGE MW	36.63kDa
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Expression Region:381-713aaStorage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:App

Note: This product is for in vitro research use only