

Cross-linked C-terminal telopeptide of Collagen alpha-1(I) chain

Catalog No: #AP79878

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Package Size: #AP79878-1 50ug #AP79878-2 100ug #AP79878-3 1mg

Description

Product Name	Cross-linked C-terminal telopeptide of Collagen alpha-1(I) chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: ICTP

Note: This product is for in vitro research use only