

## Fibroblast growth factor 1

Catalog No: #AP79879



Package Size: #AP79879-1 50ug #AP79879-2 100ug #AP79879-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Fibroblast growth factor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	FGF-1,Acidic fibroblast growth factor,aFGF,Heparin-binding growth factor 1,HBGF-1,Fgf-1,Fgfa
Accession No.	Uniprot ID:P61149
Uniprot	P61149
GeneID	25317;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Fgf1

Note: This product is for in vitro research use only