

Galanin-like peptide

Catalog No: #AP79880



Package Size: #AP79880-1 50ug #AP79880-2 100ug #AP79880-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Galanin-like peptide
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Accession No.	Uniprot ID:Q9UBC7
Uniprot	Q9UBC7
GeneID	85569;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:GALP

Note: This product is for in vitro research use only