

6-phosphofructokinase type C

Catalog No: #AP79881



Package Size: #AP79881-1 50ug #AP79881-2 100ug #AP79881-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	6-phosphofructokinase type C
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	6-phosphofructokinase, platelet type, Phosphofructo-1-kinase isozyme C, PFK-C, Phosphofructokinase 1, Phosphohexokinase, PFKF
Accession No.	Uniprot ID:Q01813
Uniprot	Q01813
GeneID	5214;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: PFKP

Note: This product is for in vitro research use only