

Pancreatic secretory granule membrane major glycoprotein GP2

Catalog No: #AP79882

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP79882-1 50ug #AP79882-2 100ug #AP79882-3 1mg

Description

Product Name	Pancreatic secretory granule membrane major glycoprotein GP2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Pancreatic zymogen granule membrane protein GP-2,ZAP75
Accession No.	Uniprot ID:P55259
Uniprot	P55259
GeneID	2813;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:GP2

Note: This product is for in vitro research use only